

Available functions	CHEESE E	LIFT	EASY	UNI	MINI .04-5	MINI .04-4
Electronic regulation of cooling water	optional	optional	optional	YES	NO	NO
Mechanical regulation of water independently for each instrument	NO	optional	optional	NO	optional	optional
Electronic regulation of cooling air	optional	NO	NO	NO	NO	NO
Reversing the rotation of MM/ Endofunction for scaler	YES	YES	YES	YES	YES	YES
Automatic blowing through the instrument	YES	YES	YES	YES	NO	NO
Turn on/ turn off on LED instruments	YES	YES	YES	YES	YES	NO
Low speed regulation of commutator MM (from 60 rpm)	YES	optional	optional	optional	NO	NO
Limitation of MM torque with an autostop function	optional	NO	NO	NO	NO	NO
Additional function (Door bell ,...)	YES	NO	NO	NO	NO	NO
Programming the parameters as: time of cup filling and bowl rinsing, time of automated switch-off of the instrument light, amount of cooling water	YES	YES	YES	YES	NO	NO
Automatic start of spittoon bowl flushing after treatment	YES	YES	YES	NO	NO	NO
OP light control from dentist and assistant table	YES	YES	YES	YES	NO	NO
Switching the FARO EVA between white and orange light (with double click on the light control button or from the dentist table, or from assistant table)	YES	YES	YES	YES	NO	NO
Automatic OP light switch off when the chair is moved	YES	YES	YES	NO	NO	NO
Safety stop control	YES	YES	YES	YES	YES	YES
Timer with memory of the last set time	YES	NO	NO	NO	NO	NO
Indication for instruments treatment by Smioil	YES	YES	YES	YES	NO	NO
Programmed parameters for # users	4 users (only with touch screen option or with SK1-08 chair)	1 user	1 user	1 user	1 user	1 user